


# 3 reasons to upgrade

Do you like to stay ahead of your game? We do too. We think your equipment should be as updated as your ideas. Now is the time to trade in your CAM 200/101/100 tensiometers for modern versions - an opportunity for you to work with our top-of-the-line instruments Theta Flex and Theta Lite. We summarized the most valuable pros of our renowned optical tensiometers.



## Software packed with features


OneAttention software has been designed with you in mind. Along with a modern and intuitive user interface it gives you the possibility to analyze your measurements live. Additionally, you get automatic contact angle and surface free energy results, as well as automatic dynamic contact angle mode.

## No cleaning with automated dispensing

The automatic disposable tip dispenser ensures you don't need to worry about cleaning as droplets are created precisely using a disposable pipette tip. The droplet volume is calculated from the real time image with advanced machine vision to ensure repeatability.

## High-end imaging for reliable results

Our Theta instruments provide camera resolution up to  $1984 \times 1264$  px with a maximum of 3009 FPS, and monochromatic LED background light. All components are carefully selected to serve you with the best possible experience and results to rely on.



Visit [biolinscientific.com](http://biolinscientific.com) for more info.